

Wizard: PEG-Ion 48 Salt Screen with PEG 1K and 4K

EB-PEG-14-B

Well	Precipitant	Salt
A1	30% (w/v) PEG 1000	0.2M Ammonium acetate
A2	30% (w/v) PEG 1000	0.2M Ammonium chloride
A3	30% (w/v) PEG 1000	0.2M Ammonium fluoride
A4	30% (w/v) PEG 1000	0.2M Ammonium formate
A5	30% (w/v) PEG 1000	0.2M Ammonium iodide
A6	30% (w/v) PEG 1000	0.2M Ammonium nitrate
A7	30% (w/v) PEG 1000	0.2M Ammonium phosphate (monobasic)
A8	30% (w/v) PEG 1000	0.2M Ammonium phosphate (dibasic)
A9	30% (w/v) PEG 1000	0.2M Ammonium sulfate
A10	30% (w/v) PEG 1000	0.2M Ammonium sulfite
A11	30% (w/v) PEG 1000	0.2M Calcium acetate
A12	30% (w/v) PEG 1000	0.05M Calcium chloride
B1	30% (w/v) PEG 1000	0.2M Diammonium tartrate
B2	30% (w/v) PEG 1000	0.2M Formic acid potassium salt
B3	30% (w/v) PEG 1000	0.2M Lithium acetate
B4	30% (w/v) PEG 1000	0.2M Lithium chloride
B5	30% (w/v) PEG 1000	0.2M Lithium citrate
B6	30% (w/v) PEG 1000	0.2M Lithium nitrate
B7	30% (w/v) PEG 1000	0.2M Lithium sulfate
B8	30% (w/v) PEG 1000	0.2M Magnesium acetate
B9	30% (w/v) PEG 1000	0.2M Magnesium chloride
B10	30% (w/v) PEG 1000	0.2M Magnesium formate
B11	30% (w/v) PEG 1000	0.2M Magnesium nitrate
B12	30% (w/v) PEG 1000	0.2M Magnesium sulfate
C1	30% (w/v) PEG 1000	0.2M Potassium acetate
C2	30% (w/v) PEG 1000	0.2M Potassium chloride
C3	30% (w/v) PEG 1000	0.2M Potassium fluoride
C4	30% (w/v) PEG 1000	0.2M Potassium iodide
C5	30% (w/v) PEG 1000	0.2M Potassium nitrate
C6	30% (w/v) PEG 1000	0.05M Potassium phosphate (dibasic)
C7	30% (w/v) PEG 1000	0.2M Potassium phosphate (monobasic)
C8	30% (w/v) PEG 1000	0.2M Potassium sulfate
C9	30% (w/v) PEG 1000	0.2M Potassium thiocyanate
C10	30% (w/v) PEG 1000	0.2M Potassium/Sodium tartrate
C11	30% (w/v) PEG 1000	0.2M Sodium acetate
C12	30% (w/v) PEG 1000	0.15M Sodium chloride
D1	30% (w/v) PEG 1000	0.2M Sodium citrate
D2	30% (w/v) PEG 1000	0.2M Sodium fluoride
D3	30% (w/v) PEG 1000	0.2M Sodium formate
D4	30% (w/v) PEG 1000	0.2M Sodium iodide
D5	30% (w/v) PEG 1000	0.2M Sodium isothiocyanate
D6	30% (w/v) PEG 1000	0.2M Sodium nitrate
D7	30% (w/v) PEG 1000	0.2M Sodium phosphate (dibasic)
D8	30% (w/v) PEG 1000	0.2M Sodium phosphate (monobasic)
D9	30% (w/v) PEG 1000	0.2M Sodium sulfate
D10	30% (w/v) PEG 1000	0.2M Sodium tartrate (dibasic)
D11	30% (w/v) PEG 1000	0.2M Tripotassium citrate
D12	30% (w/v) PEG 1000	0.2M Zinc acetate

Well	Precipitant	Salt
E1	20% (w/v) PEG 4000	0.2M Ammonium acetate
E2	20% (w/v) PEG 4000	0.2M Ammonium chloride
E3	20% (w/v) PEG 4000	0.2M Ammonium fluoride
E4	20% (w/v) PEG 4000	0.2M Ammonium formate
E5	20% (w/v) PEG 4000	0.2M Ammonium iodide
E6	20% (w/v) PEG 4000	0.2M Ammonium nitrate
E7	20% (w/v) PEG 4000	0.2M Ammonium phosphate (monobasic)
E8	20% (w/v) PEG 4000	0.2M Ammonium phosphate (dibasic)
E9	20% (w/v) PEG 4000	0.2M Ammonium sulfate
E10	20% (w/v) PEG 4000	0.2M Ammonium sulfite
E11	20% (w/v) PEG 4000	0.2M Calcium acetate
E12	20% (w/v) PEG 4000	0.05M Calcium chloride
F1	20% (w/v) PEG 4000	0.2M Diammonium tartrate
F2	20% (w/v) PEG 4000	0.2M Formic acid potassium salt
F3	20% (w/v) PEG 4000	0.2M Lithium acetate
F4	20% (w/v) PEG 4000	0.2M Lithium chloride
F5	20% (w/v) PEG 4000	0.2M Lithium citrate
F6	20% (w/v) PEG 4000	0.2M Lithium nitrate
F7	20% (w/v) PEG 4000	0.2M Lithium sulfate
F8	20% (w/v) PEG 4000	0.2M Magnesium acetate
F9	20% (w/v) PEG 4000	0.05M Magnesium chloride
F10	20% (w/v) PEG 4000	0.2M Magnesium formate
F11	20% (w/v) PEG 4000	0.2M Magnesium nitrate
F12	20% (w/v) PEG 4000	0.2M Magnesium sulfate
G1	20% (w/v) PEG 4000	0.2M Potassium acetate
G2	20% (w/v) PEG 4000	0.1M Potassium chloride
G3	20% (w/v) PEG 4000	0.2M Potassium fluoride
G4	20% (w/v) PEG 4000	0.2M Potassium iodide
G5	20% (w/v) PEG 4000	0.2M Potassium nitrate
G6	20% (w/v) PEG 4000	0.2M Potassium phosphate (dibasic)
G7	20% (w/v) PEG 4000	0.1M Potassium phosphate (monobasic)
G8	20% (w/v) PEG 4000	0.1M Potassium sulfate
G9	20% (w/v) PEG 4000	0.2M Potassium thiocyanate
G10	20% (w/v) PEG 4000	0.2M Potassium/Sodium tartrate
G11	20% (w/v) PEG 4000	0.1M Sodium acetate
G12	20% (w/v) PEG 4000	0.2M Sodium chloride
H1	20% (w/v) PEG 4000	0.05M Sodium citrate
H2	20% (w/v) PEG 4000	0.2M Sodium fluoride
H3	20% (w/v) PEG 4000	0.2M Sodium formate
H4	20% (w/v) PEG 4000	0.2M Sodium iodide
H5	20% (w/v) PEG 4000	0.2M Sodium isothiocyanate
H6	20% (w/v) PEG 4000	0.2M Sodium nitrate
H7	20% (w/v) PEG 4000	0.2M Sodium phosphate (dibasic)
H8	20% (w/v) PEG 4000	0.05M Sodium phosphate (monobasic)
H9	20% (w/v) PEG 4000	0.2M Sodium sulfate
H10	20% (w/v) PEG 4000	0.2M Sodium tartrate (dibasic)
H11	20% (w/v) PEG 4000	0.1M Tripotassium citrate
H12	20% (w/v) PEG 4000	0.2M Zinc acetate